

 VOXX LIFE

VOXX THERAPY

K N E E S U P P O R T



Introducing the newest addition
to our VoxxTherapy™ product line:
the VoxxTherapy Knee Brace.


VOXXTHERAPY | KNEE SUPPORT WITH VOXX HPT

This is a brace designed to support the knee joint while providing ergonomic protection and comfort during physical activity. Amazing for hiking, running, golf, tennis, skiing – just about any physical activity you can think of! Voxx Therapy Knee Supports are also embedded with our proprietary Human Performance Technology™. Already, users are reporting increased strength and energy with minimized pain during sports games and physical activity.



Voxx HPT is a proprietary pattern of neuro-receptor activation points called neuro-points. With Voxx HPT, users of the product realize increased overall muscular strength, better balance, and improved flexibility and stamina.

- Effective in managing knee pain caused by activities that put pressure on your knees.
- Supports muscular recovery by assisting in reducing inflammation and stiffness.
- Provides a heating effect for muscular recovery.
- Breathable fabrics provide ergonomic support.
- Suitable for all activities that put stress to the joints like running, basketball, soccer, football, golf, cycling, tennis, hiking, volleyball, skiing and much more.



VOXXTHERAPY KNEE SUPPORT SIZING CHART

Size (Inches)	Length	A		B	
		LOWER LEG MEASUREMENT		UPPER LEG MEASUREMENT	
		Relaxed	Extended	Relaxed	Extended
Small	10.6"	10"	17.5"	11.5"	20"
Medium	10.6"	11"	19"	12"	20.5"
Large	10.6"	12"	20"	13"	22"
X- Large	10.6"	12.5"	21"	14"	24"
2X - Large	10.6"	14"	23"	15"	26"

MEASURING GUIDELINES

- In order to take the circumference measurement for "A" – "Lower Leg", start at the middle of the knee and measure 5 inches down the leg and then measure the circumference. This is the lower cuff of the Knee Support.
- In order to take the circumference measurement for "B" – "Upper Leg", start at the middle of the knee and measure 5 inches up the leg and then measure the circumference. This is the upper cuff of the Knee Support.
- RELAXED refers to the Knee Support's measurement when it is not being stretched or the minimal size of the Knee Support's circumference.
- EXTENDED refers to the Knee Support's measurement when it is stretched or the maximum size of the Knee Support's circumference.
- Make sure to measure next to bare skin.
- Relaxed and Extended measurements are the range for each size of the Knee Support, make sure you fit comfortably within the range. The Knee Support needs to be a snug fit but you want to make sure it is comfortable.

Fiber Content: 75% nylon, 15% spandex, 10% elastic

WARNING: While this product may be used for injury prevention, no product can prevent all injuries. And no guarantee is made either expressed or implied that any injury will be prevented by this product.

Manufactured in a    certified facility.



VOXXTHERAPY, VOXX HPT and VOXXLIFE are registered trademarks of VoxxSports Inc.